




PEN NEEDLES

STERILE SINGLE-USE PEN INJECTOR NEEDLES
INTENDED FOR THE DELIVERY OF INSULIN





Sure and simple

PEN NEEDLES

STERILE SINGLE-USE PEN INJECTOR NEEDLES INTENDED FOR THE DELIVERY OF INSULIN

- Thin wall technology
- Double silicone coating
- Ultra sharp tip
- Universal fit

CE 0197
Class IIa

PRODUCT COMPLIANCE

- CE 0197, Class IIa Medical Devices, MDD 93/42/EEC

— ISO 9692

— ISO 11608-2

— ISO 13485

MANUFACTURING DETAILS

Legal manufacturer: Promised Medical Devices Inc., CANADA

STERILIZATION AND SHELF LIFE

Sterilization

EtO (Ethylene oxide)



Single-use only




Maximum Shelf Life


58 months



MATERIALS USED



Latex-free



PRODUCT RANGE OVERVIEW

Reference code	Length	Diameter	Color
4S-IPN-0671-31-4	4 mm	31 G (0.25 mm)	Light Green
4S-IPN-0671-32-4	4 mm	32 G (0.23 mm)	Light Green
4S-IPN-0671-31-5	5 mm	31 G (0.25 mm)	Purple
4S-IPN-0671-32-5	5 mm	32 G (0.23 mm)	Purple
4S-IPN-0671-32-6	6 mm	32 G (0.23 mm)	Blue
4S-IPN-0671-31-6	6 mm	31 G (0.25 mm)	Blue
4S-IPN-0671-31-8	8 mm	31 G (0.25 mm)	Orange



PACKAGING DETAILS

Quantity per individual box	100 units
Quantity per outer box	50 individual boxes (5000 units total)
Number of outer boxes per Euro-pallet	30

	Length	Width	Height	Weight	Quantity
Individual box	130 cm	70 cm	70 cm	0.115 kg	An individual box contains 100 units
Outer box	37.5 cm	27.5 cm	37.5 cm	6.3 kg	An outer box contains 50 individual boxes (5000 units)

Languages










5 languages on IFU, "CAUTION", and all levels of packaging:

English (EN), French (FR), German (DE), Spanish (ES), Italian (IT)

Storage conditions

Store upright in a clean, dry place. Keep away from sunlight, rain, and moisture. Do not exceed temperatures of 40°C or moisture over 80%.

Label details

	Individual box	Outer box
Product name	✓	✓
Product size	✓	✓
Quantity (pieces)	✓	✓
REF	✓	✓
LOT	✓	✓
	✓	✓
STERILE	✓	✓
CE	✓	✓
	✓	—
	✓	✓
EC REP	✓	✓
	✓	✓
	—	✓
	—	✓
	—	✓
	✓	—
	✓	—
	✓	—

For more information on products or how Nipro can help you switch patients to the 4SURE portfolio, please contact Nipro Customer Service at **0800 08 588 08**.

